



Cleaning Tool


For automated cleaning



For manual cleaning



For drying



Cleaning Tube




order code
PD-CTTU-606501
Cleaning Tube 6.0/L650
Supplied with Cleaning Holder PD-CTH



order code
PD-CTTU-600504
Cleaning Tube 6.0/L50
Supplied with Cleaning Holder PD-CTH



order code
PD-CTTU-150252
Cleaning Tube 1.5/L25
Supplied with Cleaning Adaptor IN PD-CTAD-IN



order code
PD-EZ-M
EZ Spray Adaptor
• PANA Spray Plus is not included



order code
Z182100
PANA SPRAY Plus


Primado2
Standard Attachment 300

Standard Attachment 300 80,000min⁻¹

Standard Attachment 300 is designed with new features to deliver improved visibility, smooth operation, safety, and stable rotation performance.


Variable Exposure

- Adjustable exposure lengths of 13.0/15.5/18.0/20.5/23.0/25.5mm by Twist Collet Operation




Easy Twist Collet Operation

- Bur mounting/removing and variable exposure




Universal Bur Length System

- One bur size fits all lengths of Attachment
- Easier management in operating room



Accessory

- Can slide onto the attachment




Non-Slip Surface

- Improved grip with surface etching



Smooth Connection

- Quick Joint allows for an easy attachment-motor connection



Compatible with Primado 2 Control Unit and Slim Motor



order code
P200-CU-230
Primado2 Control Unit for 230V



FC-73



FC-74

order code
FC-73/FC-74
Primado2 Single Foot Control
Primado2 Multi Foot Control



P200-SMH-HS

P200-SMH


order code
P200-SMH-HS
Primado2 Slim Motor/HS
· Hand Control Motor
· Speed : Max. 80,000 min⁻¹

order code
P200-SMH
Primado2 Slim Motor
· Speed : Max. 80,000 min⁻¹

Universal Bur Length System

Integration

One bur size fits all lengths of Attachment.





Simple Bur Mounting/Removing

The bur can be loaded/removed with simple Twist Collet Operation.


Variability

(by 2.5mm : 13.0/15.5/18.0/20.5/23.0/25.5mm)

Exposure length (6 steps) is controlled by operating the Twist Collet on the attachment. This allows for variable exposure without having to handle/change the bur.




Easy Twist Collet Operation




Lineup

Straight



Angled




Accessory



order code
P300-2-B
STD Attachment Beak

*Cannot be used with the Irrigation Nozzle.



Irrigation Nozzle



































































































order code
PD-IN-D
For P300-2S10/P300-2A10/
P300-2S20/P300-2A20
order code
PD-IN-2ES
For P300-2S40/P300-2A40

order code
PD-IN-2S
For P300-2S60/P300-2A60
order code
PD-IN-2M
For P300-2S80/P300-2A80



Cutting Accessories for Standard Attachment 300

Cutting Accessories												Accessory	STD Attachment			
											Head size					
order code	order code	order code	order code	order code	order code	order code	order code	order code	order code	order code	 (mm) Actual size					
			PDS-2RDUL-06								 0.6	 order code P300-2-B STD Attachment Beak	Straight	 order code P300-2S10 STD Attachment 10/Straight		
											 1.0			 order code P300-2S20 STD Attachment 20/Straight		
											 1.2			 order code P300-2S40 STD Attachment 40/Straight		
									PDS-2RUL-13		 1.3			 order code P300-2S60 STD Attachment 60/Straight		
											 1.5			 order code P300-2S80 STD Attachment 80/Straight		
											 1.8			 order code P300-2-B STD Attachment Beak	Angled	 order code P300-2A10 STD Attachment 10/Angled
											 2.0					 order code P300-2A20 STD Attachment 20/Angled
											 2.2					 order code P300-2A40 STD Attachment 40/Angled
											 2.5					 order code P300-2A60 STD Attachment 60/Angled
											 3.0					 order code P300-2A80 STD Attachment 80/Angled
											 3.5			 order code P300-2-B STD Attachment Beak	Angled	 order code P300-2A10 STD Attachment 10/Angled
											 4.0					 order code P300-2A20 STD Attachment 20/Angled
											 5.0					 order code P300-2A40 STD Attachment 40/Angled
											 6.0					 order code P300-2A60 STD Attachment 60/Angled
											 7.0	 order code P300-2A80 STD Attachment 80/Angled				
											 7.5					

*The exposure length of these burs are not adjustable.