


# QuadEdge Bur

CUTTING ACCESSORY




**Primado**  
—series—

# QuadEdge Bur

## New fluted bur with a unique design

- Aggressive** : Original design enables efficient cutting while ejecting bone
- Controllable**: Special bur edge limits unintentional kicking and chattering for smooth drilling
- Utility** : High-performance cutting is available across the entire bur head (top, sides, and back)



Universal Bur Length System

Description	Order code	Head size (mm)	Attachment	
QuadEdge Bur UL/3.0	PDS-2QEUL-30	3.0	for P300-2S10	P300-2A10
QuadEdge Bur UL/4.0	PDS-2QEUL-40	4.0	P300-2S20	P300-2A20
QuadEdge Bur UL/5.0	PDS-2QEUL-50	5.0	P300-2S40	P300-2A40
QuadEdge Bur UL/6.0	PDS-2QEUL-60	6.0	P300-2S60	P300-2A60
			P300-2S80	P300-2A80

- \* The exposure length of the QuadEdge bur is not adjustable.
- \* The QuadEdge bur cannot be used with Standard Handpiece 200.